

=High efficiency cleaning indicator=

# CliniCheck USW



CliniCheck USW is a new cleaning indicator with high detectability of the troubles in the operation of washer-disinfector(WD) , and with high economic efficiency. According to ISO 15883-5.

## Features

- **High economic efficiency**  
Superior cost performance enable daily use.
- **High detectability of the problems in the machine**  
Possible to judge by remains of test soil visually.
- **High usability**  
Easily attach to the basket by the special holder.
- **Easy filing**  
Because made of thin plastic sheet, file easily as the document.
- **Versatile**  
For both, conventional W/D machines and ultrasonic machines


## How to use

- ① Insert a CliniCheck USW into the holder



- ② Attach the holder to the inner side of basket.
- ③ After cleaning, check the remains of test soil and do evaluation.
- ④ File CliniCheck USW in a archiving file.
- ⑤ Recommend usage of CliniCheck USW for every cleaning.

## Archiving file




Attachment two places in one basket is effective for accurate trouble detention

Applications and inquiries

## Examples of performance evaluation



Original



Successful cycle

Always define the standard (position, detergent, temperature and process time) where the soil disappears in normal operation. Check the washing machine and conditions if any test soil appears after the cycle.

## Examples of faulty test and corrective actions

Examples of Soil spot in cleaning operation	General Corrective Actions
	<ul style="list-style-type: none"><li>● check and adjust dosage of detergent</li><li>● investigate cleaning time and temperature</li><li>● check and adjust loading conditions</li><li>● check water filter and jet nozzles</li><li>● check the inset position of CliniCheck USW</li><li>● check pump</li><li>● check tube connections and reservoir of detergent</li><li>● check pre-rinsing temp. or availability of pre-rinsing step(below 40°C) etc</li></ul>

## Packaging



- CliniCheck USW  
100 CliniCheck USW  
in PE bottle
- Holder (sell separately)

## Directions for use and storage

- Store at room temperature avoiding direct sunlight
- After use, close the lid of the bottle tightly to prevent the characteristic change of test soil in unused samples by absorption of moisture
- Don't touch the test soil spot when handling CliniCheck USW. It may influence the test result.

**Clinichem**

Hallitie 12, 47400 Kausala, Finland.

Email. [clinichem@clinichem.fi](mailto:clinichem@clinichem.fi), [www.clinichem.fi](http://www.clinichem.fi)

[www.clinichem.fi](http://www.clinichem.fi)

Tel. +358 5 38 181 80

Fax. +358 5 38 181 81