

VU30H



ULTRASOUND CLEANER 30 LITERS

[Download PDF](#)

- ✓ ultrasound frequency 38 kHz
- ✓ n. 6 transducers
- ✓ transducer power effective: 300 W / peak 1,2 kW
- ✓ field strength: 2,6 W/cm³
- ✓ sweep frequency variation function to improve treatment efficiency
- ✓ gradual frequency intensity increase function for efficient instrument treatment
- ✓ degass function for fast self degassing system to increase washing efficiency
- ✓ washing water heating function from 30 °C up to 80 °C (1,2 kW)
- ✓ heating and keeping the set temperature indication with red/green LED
- ✓ preset timer function 1, 2, 3, 4, 5, 10, 15, 20, 25, 30 min. and continuous (max 12h)
- ✓ built in drain system
- ✓ stainless steel external panels
- ✓ internal tank dimensions LxDxH: 505 x 300 x 200 mm
- ✓ external dimensions: LxDxH: 568 x 340 x 321 mm
- ✓ electrical connection: 1/N/PE 230 V - 50 Hz - 1,5 kW max
- ✓ basket VU30CSK not included

Smeg S.p.A. Via Leonardo Da Vinci, 4 - 42016 Guastalla (RE) Italy - P. IVA 01555030350
Phone +39 0522 8211 - fax +39 0522 821453 - service.instruments@smeg.it